





LIQUEEL®

SOFT SILICONE POLYURETHANE FOAM DRESSING

Performances






- ✓ absorbent and self-adhesive
- ✓ high absorption capacity
- ✓ vapor permeable outer layer
- ✓ soft polyurethane foam
- ✓ adjustable without tools

Amount of exudate (E)



- none /low 
- moderate 

Signs of infection (I)

Wound phase (P)

- Necrosis 
- Fibrin coating 
- Granulation 
- Epithelialization 
- Infection 

Wound depth (W)

- shallow 
- deep 


LIQUEEL® PT



LIQUEEL® PT is a sterile, absorbent, self-adhesive, soft wound dressing made of polyurethane foam with silicone coating. The wound dressing consists of a soft wound contact layer, a polyurethane foam layer with high absorption capacity and a vapor-permeable outer layer of polyurethane film.

LIQUEEL® PT helps maintain a moist wound environment in the presence of exudate to support the natural healing process and is indicated for the treatment of light to moderately exuding wounds.

LIQUEEL® PT is particularly suitable for the treatment of convex, concave wound situations and special body sites and shapes (e.g. heel wounds, wounds of the female breast, elbows), can be adapted to the wound situation of your patient without tools, has excellent adhesive properties and enables atraumatic dressing changes.



Performances

- ✓ absorbent and self-adhesive
- ✓ high absorption capacity
- ✓ vapor permeable outer layer
- ✓ soft polyurethane foam
- ✓ adjustable without tools

Indications

- ✓ Pressure ulcers
- ✓ Venous leg ulcers
- ✓ Arterial leg ulcers
- ✓ Non-infected diabetic Foot ulcers
- ✓ 1st and 2nd degree burns
- ✓ Scalds
- ✓ Abrasions, lacerations and small cuts



	Size	Additional information	UNIT	PZN*	REF
S	Ø 75 mm	Ø 56,5 mm bis 75 mm	1 x 10	16789030	LQ75-PT
M	Ø 150 mm	Ø 112,5 mm bis 150 mm	1 x 10	16789047	LQ150-PT
L	Ø 220 mm	Ø 165 mm bis 220 mm	1 x 10	16789053	LQ220-PT
XL	-	-	-	-	-

*The Pharmazentralnummer (PZN) is a Germany-wide identification key for medicines and other pharmacy products.






Performances

- ✓ absorbent and self-adhesive
- ✓ high absorption capacity
- ✓ vapor permeable outer layer
- ✓ soft polyurethane foam
- ✓ adjustable without tools



Amount of exudate (E)

- none /low 
- moderate 

Signs of infection (I)

- Necrosis 
- Fibrin coating 
- Granulation 
- Epithelialization 
- Infection 

Wound depth (W)

- shallow 
- deep 

LIQUEEL® GD



LIQUEEL® GD is a sterile, absorbent, self-adhesive, soft wound dressing made of polyurethane foam with silicone coating. The wound dressing consists of a soft wound contact layer, a polyurethane foam layer with high absorption capacity and a vapor-permeable outer layer of polyurethane film.

LIQUEEL® GD helps maintain a moist wound environment in the presence of exudate to support the natural healing process and is indicated for the treatment of light to moderately exuding wounds.

LIQUEEL® GD convinces by flexible application and optimal positioning on all body sites, by means of tool-free adaptation to the wound situation, excellent adhesive properties and atraumatic dressing changes.



Performances

- ✓ absorbent and self-adhesive
- ✓ high absorption capacity
- ✓ vapor permeable outer layer
- ✓ soft polyurethane foam
- ✓ adjustable without tools

Indications


- ✓ Pressure ulcers
- ✓ Venous leg ulcers
- ✓ Arterial leg ulcers
- ✓ Non-infected diabetic Foot ulcers
- ✓ 1st and 2nd degree burns
- ✓ Scalds
- ✓ Abrasions, lacerations and small cuts

	Size	Additional information	UNIT	PZN*	REF
S	5 cm x 6 cm	small	1 x 10	16789082	LQ56-GD
M	9 cm x 9 cm	medium	1 x 10	16789076	LQ99-GD
-	-	-	-	-	-
-	-	-	-	-	-

*The Pharmazentralnummer (PZN) is a Germany-wide identification key for medicines and other pharmacy products.

The logo for LF MED, featuring the letters 'LF' in a bold, sans-serif font, followed by a stylized circular symbol containing a curved line, and then the word 'MED' in a similar bold, sans-serif font.

LF MED

A large graphic element at the bottom of the page. It features a vibrant pink water splash with numerous bubbles, overlaid on a semi-transparent pie chart with several segments. The background of the entire page consists of fine, wavy, light-colored lines that create a sense of motion and depth.

*L+F Medizinprodukte GmbH
Birkenweiher Straße 4
63505 Langenselbold /Germany*

Phone: +49 6184 99 53 357

Fax: +49 6184 99 53 358

E-Mail: info@lf-med.de

www.lf-medizinprodukte.de | www.therakit.de